

# Microbiological safety cabinets - Protective functions and their limits

Thomas Hinrichs ▪ Berner International GmbH ▪ Elmshorn Germany

**Key words:** Microbiological safety cabinet, personal protection, product protection, DIN EN 12469, NSF 49, performance envelope testing, downflow, inflow

## Introduction

The effect of airflow (s. fig. 1) on the protective functions of safety cabinets will be described. [1] The protective functions in terms of personal and product protection have been tested using the microbiological test method in accordance with DIN EN 12469 [2]. The optimum performance in the form of a „Performance Envelope“ has been determined for a **BERNER FlowSafe®** safety cabinet. The variation in airflow conditions for safety cabinets reduces the performance and can lead to the loss of protective functions. Comparisons show that every safety cabinet, due to design differences, has its own PE and therefore different performance limits.

## Materials and methods

Table 1 presents the results relating to the airflow conditions and microbiological tests of personal and product protection in accordance with DIN EN 12469. Initially the Safety Cabinet (SC) was set to the specific operating point of 0,35m/s for a low turbulence downflow and to 0,44m/s for air inflow. Subsequently the flow conditions were changed gradually by 0,05 m/s in accordance with NSF 49 [3], until one or more of the limit values for personal and/or product protection had been exceeded, or if the SC's limits in terms of flow dynamics were reached, arising from design-dependent reasons.




Fig. 1: Interaction inflow and downflow within SC

## Results

**Table 1: Detection reg. the performance limits of safety cabinets against airflow variation**

**Test subject:** Microbiological safety cabinet class 2, **BERNER FlowSafe®** B-[MaxPro]<sup>2</sup>-130,  $h_A = 180$  mm; **Test specifications:** DIN EN 12469 (09.2000); DIN 12980 (06.2005), NSF/ANSI 49 (01.2007); **Location:** BERNER R&D laboratory, Elmshorn, Germany.

**Test Materials:** Nebulizer: collision nebulizer, model CN 31,  $\bar{V}_n = 0,51 \pm 0,05$  m/s,  $\Delta p = 1,4$  bar; Cylinder: stainless steel;  $d = 63$  mm; Spore suspension: *B. subtilis* var. *niger*, ATCC 9372; Agar plates: trypticase soy agar,  $d = 90$  and  $150$  mm; Thermal anemometer: Testo AG, model 454; Air capture hood: TSI GmbH model Accubalance®.

**Test Materials personal protection:**  $N_n \geq 5 \cdot 8 \times 10^8$  cfu;  $N_c \geq 300$  cfu; Six impinger: model AGI 30,  $\dot{V}_i = 12,5$  l/min,  $N_i \leq 10$  cfu; Two slit samplers: model FH 5,  $\dot{V}_s = 30$  l/min,  $N_s \leq 5$  cfu.

**Test Materials product protection:**  $N_n \geq 5 \cdot 8 \times 10^6$  cfu;  $N_c \geq 300$  cfu; 25 agar plates:  $N_a \leq 5$  cfu

Tested protective function		Personal protection		Product protection		
Replicates each single airflow combination		5		3		
Airflow combinations		42		42		
Total number of tests		210		126		
Variable		Airflow		Protective function		
No.	Downflow	Inflow	$\bar{V}_D$ in [m/s]	$\bar{V}_I$ in [m/s]	Personal	Product
0	Operating point		0,35	0,44	Yes	Yes
1	Reduced	Reduced	0,11 - 0,32	0,19 - 0,40	Yes	Yes
2	Reduced	Reduced	$\leq 0,08$	$\leq 0,18$	No	Yes
3	Reduced	Raised	0,21 - 0,33	0,43 - 0,86	Yes	Yes
4	Reduced	Constant	0,08 - 0,27	0,83	Yes	No
5	Raised	Reduced	0,40 - 0,63	0,14 - 0,40	Yes	Yes
6	Constant	Reduced	0,64	$\leq 0,13$	No	Yes
7	Raised	Raised	0,37 - 0,50	0,44 - 0,89	Yes	Yes


**Glossary:**  $h_A$ : Height of the front aperture in [mm];  $\bar{V}_n$ : Average discharge air flow velocity of the nebulizer in [m/s]; ATCC: American Type Culture Collection, Rockville, MD, USA;  $N_n$ : Minimum number of spores delivered by the nebulizer per test in [cfu]; cfu: Colony forming unit;  $N_c$ : Minimum number of spores on the control plate per test in [cfu];  $N_i$ : Limit of spores collected by six impinger per test in [cfu];  $N_s$ : Limit of spores collected by two slit sampler's per test in [cfu];  $N_a$ : Limit of spores collected with all agar plates per test in [cfu];  $\dot{V}_i$ : Volume flow of each impinger in [l/min];  $\dot{V}_s$ : Volume flow of each slit sampler in [l/min];  $\bar{V}_D$ : Average downflow velocity in the work room in [m/s];  $\bar{V}_I$ : Average inflow velocity in the work aperture in [m/s]

## Discussion


The performance limits of the protective functions are very important for determining the optimum specific operating point for every model, to guarantee the highest possible safety functions. A reduction of the airflow volume inside the SC to minimise vibration, noise level or extension of filter life results in the specific operating point being very close to the performance limit of the protective functions [4, 5]. If the airflow volume is too low, then personal as well as product safety can no longer be guaranteed (fig. 2). If the air inflow is increased and at the same time the downflow is greatly reduced, then particles from the laboratory will pass into the workspace (s. fig.1, right). If downflow is significantly higher than the inflow, then particles from the workspace will enter the laboratory (s. fig. 1, left).

In addition the data of the **BERNER FlowSafe®** SC was compared with the data from other manufacturers [4, 6] (fig. 3). The performance limits of the safety functions varied strongly between the six compared models. It shows very clearly that every model has a unique „Performance Envelope“ which is dependent on the geometry, the work aperture height, shape and size of the work surface and positioning of the specific working point.

**Fig. 2: Results of the microbiol. test of personal and product protection based on inflow and downflow.**



**Fig. 3: Comparison of Performance Envelope of different safety cabinets.**



## Summary

The performance limits of the protective functions were determined using the microbiological test method under extreme airflow velocities for a **BERNER FlowSafe®** SC. For two SC, for example, the air flow velocities could be set the same, but due to differences in their design the level of protection would not be equal. The protective functions are dependent on a multitude of design-dependent factors [6]. Summarising it could be said that the **BERNER FlowSafe®** SC, in comparison to other SC, has the widest „Performance Envelope“ and therefore guarantees the highest level of safety functions.

## Bibliography

- [1] Hinrichs, T.; Safety cabinets in clean rooms: A fundamental protective device; labor&more; Succidia Verlag; 36-37; Darmstadt; 01.2009.
- [2] DIN EN 12469; Biotechnology – Performance criteria for microbiological safety cabinets; Beuth Verlag GmbH; Berlin; 09.2000.
- [3] NSF/ANSI 49; Class II (laminar flow) Biosafety cabinet; NSF International, Ann Arbor, Michigan, USA; 01.2007.
- [4] Jones, R.J., Stuart, D.G., Eagleson, D., Greenier, J., Eagleson, J.M.; The effects of changing intake and supply air flow on biological safety cabinet performance; Appl. Occ. Env. Hyg. 1990, 5 (6), 370-376.
- [5] Kruse, R. H., Barbeito, M. S.; History of the American Biological Safety Association Part III: Safety Conferences; 1978-1987; <http://www.absa.org/abohist3.html>.
- [6] Esco Biotechnology Equipment Division; Labculture Class II Type A2 Biohazard Safety Cabinet; Esco Micro Pte. Ltd.; 2005; <http://escoglobal.com/biotech/PDF/LA216lowres.pdf>

**BERNER**

the safety system