D

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GB

User Information
Cytostatic Protective Gown and Sleeve Cover

1. Read this user information carefully!
All users should be adequately trained in the use of PPE¹⁾ annually.
¹⁾: Personal Protective Equipment

2. Area of Application and Features
Maximum protection and comfort: type tested and certified as complex PPE (category III); optimal operator and product protection (sterile version); impervious to liquids on the arms and/or front, which are coated; pleasant and comfortable to wear; material is low linting / low particle generating; knitted or elasticated cuffs at the sleeve ends; sterile and unsterile version. **Protective gown:** raised neckline; breathable back; practical velcro fastening in the neck area. **Area of application:** protective gown and protective sleeve cover for handling CMR¹⁾ drugs (e.g. cytostatic and virustatic agents). **Protective capacity:** protection from all CMR drugs or chemicals is not guaranteed. **Protective barrier:** coating that is impervious to liquids on the arms and front. Always wear with the coated side outwards. Wear sleeve cover with the seam pointing downwards. **Before use:** check for damage. Do not use damaged gown and sleeve covers.

Disposal: waste requiring supervision (waste code: 18 01 04 in accordance with 2000/532/EC); in the event of heavy contamination, waste requiring special supervision (waste code: 18 01 08* in accordance with 2000/532/EC); collect and dispose of waste separately.
¹⁾: Carcinogenic, Mutagenic, toxic to Reproduction

3. Types

3.1. Protective gown

Size	S	M	L
Chest (cm)	76-92	92-108	108-124
Height (cm)	152-164	164-176	176-194

3.1.1. Blue gown with knitted cuffs

Item no.	unsterile	pu	15	item	BI-6700	BI-6800	BI-6900
	sterile		10		BI-6701	BI-6801	BI-6901

3.1.2. Light blue gown with elasticated cuffs

Item no.	unsterile	pu	15	item	BI-6500	BI-6550
	sterile		10		BI-6600	BI-6650

3.2. Protective sleeve cover

Size	Universal (length app. 520 mm)
Item no.	BI-6000
	BI-6001

3.2.1. Blue sleeve covers with knitted cuff

Item no.	unsterile	pu	50	pair	BI-6000
	sterile		40		BI-6001

3.2.2. Light blue sleeve covers with elasticated cuff

Item no.	unsterile	pu	50	pair	BI-6200
	sterile		30		BI-6300

4. Sterilisation
Fumigation with ethylene oxide.

5. Material

Material	Spun polypropylene
Material weight	42 g/m ²
Liquid-tight coating	Polyethylene (thickness approx. 25 µm)
Total weight of gown / pair of sleeve covers	120-141 g / 20-36 g

5.1. Mechanical Properties of Material

Requirement	Performance Level
Abrasion resistance (1-6) acc. EN 530	2
Puncture resistance (1-5) acc. prEN 863	1
Seam resistance (1-5) acc. ISO 5082	2
Tear propagation resistance (1-5) acc. ISO 4674	Longitudinal 2 Transverse 3

5.2. Protection from Chemical Hazard; permeation¹⁾ acc. EN 374 part 3

Chemical	Breakthrough Time ²⁾ [min]	Performance Class(1-6) ³⁾
Carmustine	> 140	4
Amsacrin	> 480	6
Cisplatin	> 480	6
Cyclophosphamide	> 480	6
Doxorubicin	> 480	6
5-Fluorouracil	> 480	6
Methotrexate	> 480	6
Paclitaxel	> 480	6
Thiotepa	> 480	6
Vincristine	> 480	6
NaOH 30%	> 480	6

¹⁾: Movement of a chemical through a material on a molecular level. ²⁾: At a permeation rate of 1 µg/min·cm². ³⁾: The performance class does not reflect the actual duration of protection at the workstation!

6. Proper use
Protective gown/cuff replacement interval: daily, i.e. use for a maximum of 8 hours; immediately in the event of visible contamination. Do not reuse!

7. Care instructions
Do not wash, do not iron, do not tumble dry, do not dry clean.

8. CE Mark
CE mark for complex PPE in category III in accordance with PPE Directive 89/686/EEC, based on DIN EN 467 (01.99); DIN EN 374 part 3 (04.94); DIN EN 340 (08.93). Type test and control measures by intermediary notified body "0299".

9. Notified Body "0299"
Technical Committee for Personal Protective Equipment, testing and certification body of BG-PRÜFZERT Centre for Safety Technology, Zwengenberger Str. 68, D-42781 Haan, Germany.

10. Quality management system
Tested and certified by TÜV Product Service GmbH in accordance with DIN EN ISO 9001:2000.

11. Storage and transport
Dark (protect from direct UV light and sunlight); cool (+5 to +40°C); dry; no contact with pointed and/or sharp objects.

12. Date of manufacture
1.+2. number = month, 3.-6. number = year.

13. Shelf life
Unsterile version: 5 years from the date of manufacture.
Sterile version: 4 years from the date of sterilisation.

14. LOT number
LOT number unsterile / sterile version: 6 / 4 numbers, numbered consecutively.

15. Copyright and industrial property rights
The **Z*** product group sign is a registered and protected trademark. Infringement of industrial property rights will be prosecuted. BERNER INTERNATIONAL GMBH assumes no warranty or liability for protective gown and sleeve cover that are employed beyond the scope of their intended use. The German version of the user information is binding for translations into foreign languages. The content of the user information can be changed at any time and without notice. This user information is the property of BERNER INTERNATIONAL GMBH and may not be forwarded to third parties, copied, photocopied or used for promotional purposes, either wholly or in part, without written consent. Product specifications are subject to change. © BERNER INTERNATIONAL GMBH

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