



ISYSOFT Cleaning Cloths (sterile / non sterile)
ISYSOFT Cleaning Cloths latex-free (sterile / non sterile)



Area of Application and Features

- **Area of application:** Cloths for cleaning, especially suitable for cleaning surfaces contaminated with CMR¹⁾ drugs
- **Properties:** Extremely absorbent cloths, folded several times, suitable for the daily cleaning routine and for absorbing spilled liquids. The cloths can be used dry or moistened/wet. ISYSOFT-cleaning cloths are low in lint and are highly effective and are highly effective at cleaning. For single use only!
- **Recommended use:** For the cleaning of cytostatic contaminations we recommend, based on a Study by the German Institute for Energy and Environmental Technology (IUTA e.V.):

Cleaning in 2 steps:

first using a) an alkaline cleaning agent > pH 11 and b) Isopropanol (70%)

In a study on the subject of cleaning cloths carried out by the Institute for Energy and Environmental Technology, highly contaminated surfaces (with different cytostatics) were cleaned using ISYSOFT cleaning cloths, which reached a degree of cleanliness of 99.95%.

Wipes are moistened with the cleaning solution, avoiding spray. In case of highly contaminated surfaces (e.g. after spills, etc. the 2-step process should be repeated 3 times.

Disposal: If cloths have been used to clean contaminated areas (CMR-drugs): waste requiring supervision (waste code: 18 01 04 in accordance with 2000/532/EC); in case of high contamination, waste requiring special supervision (waste code: 18 01 08* in accordance with 2000/532/EC); collect and dispose of waste separately.

¹⁾: Carcinogenic, mutagenic, reprotoxic


ISYSOFT (non sterile version)

Dimensions	400 x 300 mm
Item-no.	121500
Colour	natural white
Packaging unit	1 PU = 100 pieces

Material Properties

Material	airlaid paper
Material composition	cellulose and binding agent made of latex
Material weight	80 g/m ²
Folding	twice
Absorbable liquid volume¹⁾	≤ 65 ml

¹⁾: Internal tests.

ISYSOFT (sterile version)	
Dimensions	400 x 300 mm
Item-no.	121501
Colour	natural white
Packaging unit	1 PU = 50x 2 pieces (ideal for 2-step cleaning)
Material Properties	
Material	airlaid paper
Material composition	cellulose and binding agent made of latex
Material weight	80 g/m ²
Folding	three times
Absorbable liquid volume¹⁾	≤ 65 ml
¹⁾ : Internal tests.	
 Sterilisation	
Procedure	fumigation with ethylene oxide

ISYSOFT latex-free (non-sterile version)	
Dimensions	400 x 300 mm
Item-no.	121505
Colour	white
Packaging unit	1 PU = 100 pieces
Material Properties	
Material	spunlace non-woven fabric, latex-free
Folding	twice
Absorbable liquid volume¹⁾	≤ 78 ml
¹⁾ : Internal tests.	

ISYSOFT latex-free (sterile version)	
Dimensions	400 x 300 mm
Item-no.	121506
Colour	white
Packaging unit	1 PU = 50x 2 pieces (ideal for 2-step cleaning)
Material Properties	
Material	spunlace non-woven fabric, latex-free
Material weight	50-55 g/m ²
Folding	three times
Absorbable liquid volume¹⁾	≤ 78 ml
¹⁾ : Internal tests.	
Sterilisation	
	
Procedure	fumigation with ethylene oxide
Quality Management System	
	
Our quality management system is tested and certified by TÜV Product Service GmbH (a certification body accredited by the German Accreditation Council,) in accordance with DIN EN ISO 9001:2008. Regular audits and production site inspections guarantee the quality of our products.	
Storage and Transport Conditions	
	
<ul style="list-style-type: none"> • Dark (protect from direct UV light and sunlight) • Cool (-5 to +40°C) • Dry • No contact with pointed and/or sharp objects 	
Shelf Life	
	
<ul style="list-style-type: none"> • 3 years from the date of manufacture 	
Barcode-system for storage	
<ul style="list-style-type: none"> • On all packaging units • Item-nos. are readable as EAN-13 code • Manufacturing and use-by date as well as lot number (LOT) is readable as code 39 	